Lampa foto Refine MaxCure9

Pret:

663,00 lei TVA inclusa

Descriere produs:

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Broad-spectrum
AI Intelligence 1 sec cure

MaxCure9
Curing 2mm Resin in 1s

01 Dual intensity mode


- P1 High intensity mode


Curing time
For orthodontic brackets, post and core and curing photosensitive adhesive in porcelain veneer

- P2 General intensity mode

Curing time
For Normal resin

Meet Different Clinical Needs


0 3 5 Time (S)

High intensity mode

Constant intensity: 2300-2500mw/cm²
Used for quick curing, ensuring enough curing depth, more efficient
Time setting: 1s, 3s

General mode

Constant intensity: 1000-1200mw/cm²
for most treatment scenarios like teeth filling, restoration etc.,
Time setting: 5s, 10s, 15s, 20s


02 Broad-spectrum

Curing most materials

Most of the resins on the market use CQ as the photoinitiator, but some resins use two new types of photoinitiators, diphenylphosphorus oxide or Ivocerin, and the new materials containing TPO have a light absorption peak of 405nm. If you use an ordinary narrow-spectrum light curing machine (wavelength range: 420nm-480nm), the resin is prone to incomplete curing. And we know that incomplete composite curing often leads to secondary pulpitis, apical inflammation, and tooth hypersensitivity.

Broad-spectrum

Curing most materials

MaxCure 9 broad-spectrum light curing, newly added the 385-420nm band, the wavelength range is 385nm-515nm. It can better cure resin materials containing TPO and other emerging photoinitiators, and is suitable for ideal curing of all resin materials on the market


After a full upgrade, MaxCure 9 changes from a single blue light to a blue-violet light, which can truly cure all light-curing materials

03 Super large battery capacity

Battery capacity Comparison among four Brands


lithium battery- International CB certification
MaxCure9 Battery capacity is 85% higher than normal Curing Lights

04 AI Intelligence system

Advanced technology

A CONSTANT CURRENT DRIVER

Constant light intensity
Constant light intensity output, stabilizing the working power

D TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com. : I/09/137/09 Cod fiscal: RO 117 58 273

LIGHT GUIDE AUTOMATIC DETECTION

The main unit will detect the light guide and stop the "Er" flashing if the guide is connected well, or the main unit will close down after 10s flashing.

Well-installed during the period of the error alarm will stop the error alarm and the lamp will light up instantly when the lamp holder is installed.

OVERHEATING PROTECTION

Main unit will detect the temperature. The overheating protection program will run when the curing light reaches a certain temperature after long time use. It'll keep flashing until 200s later that it recovers from high-temperature.

BATTERY PERFORMANCE DETECTION

When detect an inconstant power output from battery, the main unit will stop working to warn with flashing "Er"

OVERVOLTAGE PROTECTION

With overvoltage adapter, MaxCure 9 is under protection and can't be charged. The original adapter is recommended.

Charging

fully charged

Blue light

white light

MaxCure9
Curing 2mm Resin in 1s

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273