

Aspirator Extraoral VacStation



Pret:

3.059,57 lei TVA inclusa

Descriere produs:

Detalii suplimentare: 0722 350 863 / 0729 994 382

Inscrieri pentru precomanda: contactand agentul dvs. Tehnical Dent sau prin email la comenzi@technicaldent.ro.

[PREZENTARE EXHAUSTIVA Q&A VACSTATION](#)

VacStation este un sistem extraoral de vacuum al aerosolilor care are rolul de a preveni raspandirea sangelui, virusilor, prafului si a altor substante in cabinetul dentar, in urma proceselor de frezare cu turbina, laser dentar, electrochirurgie, etc.

Prezinta 3 sisteme de filtrare foarte eficienta, precum si un sistem de sterilizare cu raze ultraviolete.

Nu prezinta un zgomot mare si are o putere de absorbtie foarte mare a aerosolilor datorata motorului cu inductie, fara perii, cu o durata de viata foarte mare.




Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB

0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

LISTA DE PRETURI SI PERIOADE DE SCHIMB FILTRE SI CONSUMABILE

D TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehncaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com : I/09/137/99 Cod fiscal: RO 117 58 273





HEPA filtration system
H13 level, trap virus and
germs $\geq 0.3\mu\text{m}$ with 99.97% efficiency

High-fiber moisture intake cotton filter
aerosol capture and moisture drying

HEPA filter
catch dust and spit with 99% efficiency

Efficient compound filter
(include another HEPA H13 filter)
with the function of adsorption, sterilization
and filtration



Brushless Motor
higher power ,longer life, lower noise

UV light
intensity 1210uw/ cm2 , wavelength 280nm,
killing virus and bacteria with high efficiency




Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273


VacStation

Eighteenth

The 1st high-fiber filter



The 2nd HEPA 13 filter



The 3rd Efficient compound filter system







Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: I/09/137/99 Cod fiscal: RO 117 58 273



VacStation is an extraoral aerosol vacuum device to prevent blood, virus, dust and other substance from spreading in dental treatment like high-speed drilling, dental laser and electro surge procedures, etc.



D TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068.0043.7192.0001 Nr. Inreg. Reg. Com : I/09/137/99 Cod fiscal: RO 117 58 273



Model	VacStation
Voltage	AC220V/50Hz
Power	250W Max
Maximum Suction Power	3000Pa (10 levels adjustment)
Maximum Volume	3.7m ³ /Min
Noise Decibel	62db
Weight	21kg
Dimensions	70cm*34cm*65cm



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Test putere absorbtie VacStation:

Unboxing Vacstation Eighteeth

Demonstratie folosire VacStation

Metoda de asamblare VacStation