

Sealer bioceramic Well-Root ST

Pret:

428,40 lei TVA inclusa



Descriere produs:

Foarte potrivita atat pentru coafaje directe cat si indirecte!

Well-Root ST™ este o pastă bioceramică injectabilă, preamestecată pentru a fi ușor de folosit, gata de utilizare, proiectată pentru obturarea permanentă a canalelor radiculare.

Well-Root ST™ are o compoziție de silicat de calciu, care necesită prezența umidității pentru a se întări.

Well-Root ST™ nu se contractă în timpul prizei și demonstrează proprietăți fizice excelente.

- Aderenta excelenta la gutaperca si dentina
- Con tracție zero în timpul întăririi
- Generează hidroxiapatită

Prezentare

- 1 x seringă de 2 g și 15 vârfuri de unică folosință (20 GA).

TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB

0068 0043 7192 0001 Nr. Inreg. Reg. Com. : I/09/137/09 Cod fiscal: RO 117 58 273



• **Tight Sealing**

- Excellent bonding to gutta-percha and dentin
- Zero shrinkage
- Hydroxyapatite generation
- Hydrophilic
- Extremely small particle size

► **Right Image**


Optical microscope and SEM image of GP/Well-Root ST /Dentin interface

: **No gap at interface**




Optical microscope image of microleakage test (dye penetration test)
; **No dye penetration**

• **Radiopacity**



Well-Root ST uses zirconium oxide for radiopaque agent.

• **High pH**



- Calcium hydroxide is generated during setting
- pH value increases 12 and over
- Alkaline environment making it highly antibacterial