

Autoclav 16 litri clasa N



Pret:

4.375,01 lei TVA inclusa

Descriere produs:

Autoclav clasa N (fara vacuum) destinat cabinetelor stomatologice, cabinetelor medicale, saloanelor de cosmetica.

Camera de sterilizare din otel inoxidabil, dintr-o singura bucata fara suduri sau imbinari, cu o durata de viata foarte indelungata!

Suport tavite si tavite integral fabricate din otel inoxidabil. Doua modalitati (directii) de plasate a suportului de tavite.

Modalitate foarte usoara de curatare a garniturii de usa.

Display digital al temperaturii cu o deviere a temperaturii setate de doar 0.5 C.

Putere de incalzire de 1800W, foarte puternica, pentru a face intreg procesul de sterilizare foarte rapid!

[SFATURI SI INSTRUCIUNI REFERITOARE LA AUTOCLAVUL CLASA N](#)

D TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@technicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Date tehnice:

- Alimentare la priza 220V
- Putere 1800W
- Capacitate camera de sterilizare: 16 litri
- Timpi de sterilizare (excluzand preincalzirea): 121C (20 minute), 134C (6 minute)
- Dimensiune camera sterilizare: diametru 230 mm x 360 mm adancime
- Dimensiune produs: 440 x 575 x 370 mm

SEA-16L-B-LED

**Within Touch
And Endlessly Clean**





PRODUCT FEATURES

- SS304 chamber, one piece, without welding
- SS304 tray and rack. The rack could be in two directions
- Gasket is easy to be cleaned
- Temperature digital display
- Temperature deviation is less than 0.5°C
- 1800W heating elements, speed up the heating, more efficient

SEA-16L-B-LED

Basic Technical Data:

Voltage: 220V or 120V

Frequency: 50Hz/60Hz

Power: 1 800W

Chamber Capacity: 16L

European Class N

Sterilization temperature:

121 °C/20min, 134 °C/6min

Chamber Size: Ø230x360mm

Product Size: 440x575x370 mm

