

Ansa detartraj Ultramint G1



Pret:

25,89 lei TVA inclusa

Descriere produs:

Ansa detartraj din material rezistent, inoxidabil, cu filet interior, compatibila cu aparatele de detartraj Ultramint si Ultramint Pro, Woodpecker, EMS

Indicatii:

Pentru indepartarea depunerilor supragingivale in toate cadranele.

Model: G1

TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehncaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com. : I/09/137/09 Cod fiscal: RO 117 58 273






EMS **G1**

Order Code **6680523**

Satelec **G1-S**

Order Code **6680617**

Function Description: All-round removal of supragingival calculus



Mode: G
Power: 1-10



Recommended
Service Life
200 Teeth



Water Flow Rate
mid&high