

Freza diamantata Diatech G379



Pret:

90,77 lei TVA inclusa

Descriere produs:

Instrumente diamantate Diatech pentru prepararea cavitatii.

Diametru: 016, 023

Granulatie:

Negru - extradura

Albastru - dura

D TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB

0068 0043 7192 0001 Nr. Inreg. Reg. Com. : I/09/137/09 Cod fiscal: RO 117 58 273

Rosu - medie

Galben - fina

Alb - ultrafina

Mod de prezentare:

Set 5 freze diamantate Diatech

G379 Egg



ISO 806...277...014



ISO 806...277...016



ISO 806...277...018



ISO 806...277...023



ISO 806...277...023



ISO Ø 1/10 mm	014	016	018	023	023
Head L in mm	3	3.5	3.5	4.2	5