

Freza diamantata Diatech G878K



Pret:

90,77 lei TVA inclusa

Descriere produs:

Instrumente diamantate Diatech pentru prepararea cavitatii.

Granulatie:

Negru - extradura

Albastru - dura

Rosu - medie

Galben - fina

Alb - ultrafina

Mod de prezentare:

Set 5 freze diamantate Diatech

G878K Tapered Torpedo



ISO 806...298...010



ISO 806...298...012



ISO 806...298...014



ISO 806...298...016



ISO 806...298...018



ISO 806...298...021



ISO Ø 1/10 mm	010	012	014	016	018	021
Head L in mm	8	8	8	8	8	8