

Camp steril 50 x 70 cm



Pret:

8,29 lei TVA inclusa

Descriere produs:

Camp steril 50 x 70 cm

Indicatii:

- Campuri sterile Steri Quick pentru acoperirea tavilor
- Protejeaza pacientul
- Acopera si protejeaza instrumentelor sterille

Beneficii:

- Fiecare camp steril este impachetat individual steril si gata de folosire

Prezentare:

- 1 x camp steril 50 x 70 cm



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273