

## HERAENIUM PW

Pret:

**1,60 lei TVA inclusa**



### Descriere produs:

**Compozi?ie:** Co 55%, Cr 24%, W 15%, Si 1%, Mn 0,8%, Fe 4%.

### Caracteristici:

- Densitate 8,2g/cm<sup>3</sup>
- Coeficient de dilatare termica 14,3μm/mK (25 – 500 grade C)
- Duritate 290 HV
- Interval de topire: 1320 grade C-1400 grade C
- Temperatura de turnare 1550 grade C