

## Gips sintetic alb clasa a III-a

Pret:

**118,14 lei TVA inclusa**



Descriere produs:

Raport ap? / pulbere	Setarea timpului	Setarea expansiunii 2h	Rezisten?a la compresiune 1h	Rezisten?? la compresiune uscat?
23 ml / 100g	10 - 13 min.	? 0,11%	40 MPa (5800 psi)	77 MPa (11165psi)

Cantitate:

5 KG



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: [office@tehnicaldent.ro](mailto:office@tehnicaldent.ro) IBAN: RO17 RNCB  
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273