



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

## Cure Box Dual cuptor foto

**Pret:**

**2.168,00 lei TVA inclusa**



### **Descriere produs:**

Dispozitiv de intarire a rasinii.

### **Caracteristici**

Lungime de undă selectabilă (400 ~ 405nm / 460 ~ 465nm)

Polimerizare mai rapidă folosind lampa LED - luminozitate de 4 ori mai puternică

Alimentare stabila

Luminozitate reglabila

### **Specificatii tehnice**

Putere: AC 100V-240V, 50/60Hz

65W



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

Interval de temperatura: 0°C~40°C

Greutate: 1,2 kg

Umiditate: 50%~80%

### **Mod de prezentare**

1 x Cure Box Bual