

## Cure Box Dual cuptor foto



**Pret:**

**2.248,53 lei TVA inclusa**

### **Descriere produs:**

Dispozitiv de intarire a rasinii.

### **Caracteristici**

Lungime de und? selectabil? (400 ~ 405nm / 460 ~ 465nm)

Polimerizare mai rapidă folosind lampa LED - luminozitate de 4 ori mai puternică

Alimentare stabila

Luminozitate reglabila

### **Specificatii tehnice**

Putere: AC 100V-240V, 50/60Hz

65W

Interval de temperatura: 0°C~40°C

Greutate: 1,2 kg

Umiditate: 50%~80%



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: [office@tehnicaldent.ro](mailto:office@tehnicaldent.ro) IBAN: RO17 RNCB  
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

### **Mod de prezentare**

1 x Cure Box Bual