

Model practica educationala silicon sutura laparoscopica



Pret:

42,53 lei TVA inclusa

Descriere produs:

Model practica educationala silicon pentru sutura laparoscopica

Prezinta un singur strat din silicon

Model cu diferite modele

Marginile sunt iesire in relief

Lungime: 13,5 cm

Latime: 10,8 cm

Grosime: 4 mm

Material: silicon hipoalergenic

Mod de prezentare:

1 x model practica educationala



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273