


Linguri de amprenta autoclavabile



Pret:

31,57 lei TVA inclusa

Descriere produs:

- Linguri de amprenta autoclavabile, pot fi incalzite la 121 grade C, se pot refolosi dupa sterilizare.
- Design unic din plastic rezistent si durabil.
- Design ergonomic, suprafetele netede asigura confort pacientului.

Prezentare

- punga 10 buc asortate



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273