

Tavita mica rasina pentru imprimanta Accufab L4D Shining 3D

Pret:

1.779,41 lei TVA inclusa



Descriere produs:

Tavita de volum mic de rasina, echipata cu o noua tehnologie a filmului, imbunatateste acuratetea si fiabilitatea tiparirii. Acest rezervor compact de rasina se mandreste cu o longevitate exceptionala, sustinand pana la 120.000 de straturi si permitand sesiuni de imprimare extinse cu numeroase printari de succes.

Dimensiune: 313.8 mm x 204. 85mm x 44mm

Volum: min 200ml - max 350ml



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273