

## Oferta Oxasil Kit silicon condensare

Pret:

**216,20 lei TVA inclusa**



### Descriere produs:

Silicon de condensare ce prezinta o prepolimerizare inainte de amestec, oferind un efect "snap-set" similar siliconilor de aditie, garantand un timp de amestec si de lucru fara prepolimerizare.

Stabilitate dimensionala 7 zile.

### Indicatii

- pre-amprentare pentru tehnica "double mix"
- amprentarea "sandwich"
- amprente functionale
- timp scurt de priza intr-orală
- reproducerea precisa a detaliilor prin Hydrocontrol



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: [office@tehnicaldent.ro](mailto:office@tehnicaldent.ro) IBAN: RO17 RNCB  
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

**Prezentare**

- Kit continand:
- oxasil 900ml
- Oxasil tube light flow 140 ml
- Oxasil activator