

D TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@technicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Primsole Light Body Fast 2x50ml



Pret:

163,19 lei TVA inclusa

Descriere produs:

Materiale de amprenta (VPS) 2 x 50ml

AVANTAJE

Rezistent: rezistenta ridicata la rupere

Reproducere inalta a detaliilor

Stabilitate dimensionala lunga

Recuperare elastica excelenta

Siguranta

Termen de valabilitate de 3 ani

Combinat cu varfuri de amestecare avansate - usor de blocat

CARACTERISTICI

Consistenta tixotropica



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273
Corp usor tip 3 - ISO 4823

Hidrocompatibil

Timp rapid in cavitate

Biocompatibil

Scanabil

Culoare: dark blue

Prezentare

Primsole Light Body Fast - 2 x 50ml, 12 varfuri