

Primsole Light Body Normal 2x50ml

Pret:

158,10 lei TVA inclusa



Descriere produs:

Materiale de amprenta (VPS) 2 x 50ml

AVANTAJE

Rezistent: rezistenta ridicata la rupere

Reproducere inalta a detaliilor

Stabilitate dimensionala lunga

Recuperare elastica excelenta

Siguranta

Termen de valabilitate de 3 ani

Combinat cu varfuri de amestecare avansate - usor de blocat



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

CARACTERISTICI

Consistentă tixotropica

Corp usor tip 3 - ISO 4823

Hidrocompatibil

Timp rapid în cavitate

Biocompatibil

Scanabil

Culoare: light blue

Prezentare

Primsole Light Body Normal - 2 x 50ml, 12 varfuri