


## Tavita Aluminiu Instrumentar Autoclav SEA-23L Runyes

**Pret:**

**171,01 lei TVA inclusa**



### Descriere produs:

### Caracteristici:

- **Material:** Aluminiu rezistent la temperaturi înalte
- **Grosime:** 2 mm
- **Greutate:** 250 grame
- **Dimensiuni:** 28 cm x 18 cm

### Compatibilitate:

- Runyes SEA-23L
- Tecno-Gaz Europa B Evo 24L
- Dr. Mayer Class B 23L