

Servetele suprafete Eurosept Xtra Refill XS



Pret:

35,01 lei TVA inclusa

Descriere produs:

Servetele dezinfectante pentru suprafete pe baza de alcool, pentru suprafete rezistente la alcool. Nu utilizati pe acrilice.

Timp de contact (raportat la lichidul dezinfectie):

bacterii - EN13727, EN1276, EN13697, EN16615: 1 min;

Candida albicans - EN13624, EN1650, EN13697, EN16615: 1 min;

TB (tuberculoza), agenti patogeni (M. terrae) - EN14348: 1 min;

micobactericid (Mycobacterium avium) - EN14348: 1 min;

virusuri invelite (HBV, HCV si HIV) - EN14476, EN16777: 1 min;

norovirusuri: 1 min (EN14476, clean), 2 min (EN16777, clean);

adenovirusuri: 1 min (EN14476, EN16777, clean);

TEHNICALDENT

Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273
polyoma SV40 (RKI/DVV dirty): 2 min;

virucid limitat plus (adenovirusuri, norovirusuri, rotavirusuri) - EN14476: 1 min (clean);

virucid limitat plus (adenovirusuri, norovirusuri, rotavirusuri) - EN16777: 2 min (clean).

Timp de contact raportat la servetelele impregnate:

bacterii si Candida albicans - EN16615 (contaminare ridicata): 5 min.

Fabricat din material textil mixt PP/VIS (12 x 20 cm).

Dimensiune: XS - 12 x 20 cm

Prezentare

Pachet 120 bucati