

Compresor dentar Energy F114 24/114



Pret:

3.681,85 lei TVA inclusa

Descriere produs:

- Putere electrica: 750W / 1 CP
- Presiune maxima de lucru: 8 bar
- Presiune normala de lucru: 5-7 +/- 0,5 bar
- Debit de aer: 114l/min
- Capacitate rezervor: 24l
- Nivel de zgomot: 70 dB
- Dimensiuni: 41 x 40 x 68 cm (latime x adancime x inaltime)



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB

0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

- Masa: 25 Kg
- 1400 rpm
- 1 cilindru
- Potrivit pentru 1 unit